

Testing Supply Substitution Strategies

This resource is intended for labs performing COVID-19 tests that are authorized. This resource includes validated supply alternatives that labs can use to continue performing testing when there is a supply issue with some components of a test.

The information in this resource is not intended to alter any already issued EUA for a COVID-19 diagnostic test nor is it intended to speak to any specific FDA regulatory requirement. Rather, the information is being provided to help address availability concerns regarding certain critical components of COVID-19 diagnostic tests during this pandemic.



[START: Go to the Main Menu](#)

Real-time RT Polymerase Chain Reaction (PCR) Component Substitution Strategies

Select a button to learn more about a topic

Specimen Collection

Learn about specimen collection, swabs, and media

Intro to PCR

Learn about RT-PCR: what it is and what it's used for

System Types

Compare use of Open Mix and Match Systems vs. Closed Systems, including instruments, supplies, and reagents

PCR Testing Process

Learn about the steps in the PCR testing process

Substitution Options

Learn about substitution options for the Open Mix and Match style CDC test

Return to the Main Menu

Additional Resources


Return to FDA FAQ

Specimen Collection

Specimen collection is the process of obtaining a sample from a patient, usually by swabbing the nose or mouth, then placing the swab in a tube that is commonly filled with liquid (media) which maintains the sample for transport to the lab.

Originally, specimen collection for SARS-CoV-2 testing required a specialized Nasopharyngeal (NP) swab and Viral Transport Media (VTM).

Options for swabs and media are below.



Explore acceptable choices for swabs

Explore acceptable choices for media

For more information please see the [FDA FAQ on this topic](#) under “Testing Supply FAQ.” It is important that the swab type be appropriate for the anatomic site on which it is used, i.e. only a Nasopharyngeal swab should be used to obtain a Nasopharyngeal specimen.”


Return to the Main Menu

Additional Resources

Return to FDA FAQ

Specimen Collection: Swabs

Choices for Swabs*			
NasoPharyngeal (NP) Swab	OroPharyngeal (OP) Swab	Mid-Turbinate (MT) Swab <ul style="list-style-type: none"> Flocked, tapered swab 	Anterior Nares (“Nasal”) Swab <ul style="list-style-type: none"> Round Foam Spun Polyester



* For more information please see the [FDA FAQ on this topic](#) under “Testing Supply FAQ.”
 It is important that the swab be appropriate for the anatomic site on which it is used, and that the swab type (e.g. polyester vs rayon) is compatible with that platform. Rayon swabs may not be compatible with all molecular testing platforms. Analytical testing should be performed to confirm compatibility with individual platforms.

Return to Specimen Collection

Return to PCR Testing Process
 “Collect Specimen from Patient”

Return to the Main Menu

Additional Resources

Return to FDA FAQ

Specimen Collection: Media

Choices for Media			
VTM <ul style="list-style-type: none"> e.g., Universal Transport Media (UTM) 	Liquid Amies-based Media <ul style="list-style-type: none"> e.g., Eswab 	Inactivating Transport Media <ul style="list-style-type: none"> e.g., Primestore MTM* 	Saline Solution <ul style="list-style-type: none"> Normal saline Phosphate-buffered saline (PBS)



As explained in FDA’s [COVID-19 Transport Media Policy](#), all transport media, including saline, should be sterile to avoid contamination of the specimen.

WARNING: Media containing guanidine thiocyanate or similar chemicals produce a potentially hazardous chemical reaction that releases cyanide gas when exposed to bleach.

Please see the full Warning statement related to this issue for more information, including a list of some media and devices where extra care should be taken because of this issue.

See **WARNING Statement**

Return to Specimen Collection

Return to PCR Testing Process
“Collect Specimen from Patient”

Return to the Main Menu

Additional Resources

Return to FDA FAQ

WARNING - Specimen Collection: Media

Many inactivating transport media (ITM) contain guanidine thiocyanate or similar chemicals making them incompatible with certain in vitro diagnostic platforms, including those that use bleach (sodium hypochlorite), and with laboratory processes that use bleach. See warning below.

WARNING: Media containing guanidine thiocyanate or similar chemicals, including Longhorn PrimeStore MTM, Zymo DNA/RNA Shield, as well as that found in the Spectrum Solutions Saliva Collection Devices, should not be used with in vitro diagnostic platforms that use bleach (sodium hypochlorite) in the testing platform or during laboratory processes.

Media containing guanidine thiocyanate or similar chemicals produce a **potentially hazardous chemical reaction that releases cyanide gas when exposed to bleach**. Molecular platforms that use bleach in a disinfecting step include the **Hologic Panther** and **Panther Fusion** systems.

Additionally, media containing guanidine thiocyanate or similar chemicals should not be used with tests that are not designed for use with such chemicals. If you have questions about compatibility with guanidine thiocyanate media and your EUA test or test platform, please contact the test or platform manufacturer for further recommendations.

**BACK to Specimen Collection:
Media**

What Is RT PCR?

RT-PCR stands for “Reverse-transcriptase polymerase chain reaction” (commonly referred to as “PCR”)

RT-PCR is a method used to detect RNA nucleic acids (a type of genetic material)

- PCR is used to detect the presence of SARS-CoV-2, the virus that causes COVID-19 disease
- Finding SARS-CoV-2 genetic material in a specimen indicates that a person is infected with the virus

Overview of System Types

Select a system type for more information

Closed System



Extraction and Amplification

Open Mix and Match

FDA



Extraction



Amplification

Return to the Main Menu

Additional Resources

Return to FDA FAQ

Released 7/23/2020

Closed System

Instruments

A single dedicated instrument performs extraction and amplification



Extraction and PCR System

Supplies and Reagents

Some materials may be specific to the instrument platform




Proprietary Reagent Kit



Collection and Media



Proprietary Reagent Cartridge



General Reagents

[Return to the Main Menu](#)

[Return to Overview of System Types](#)

[Additional Resources](#)

[Return to FDA FAQ](#)

Open Mix and Match (CDC Design)

Instruments

Extraction and amplification performed separately



Extraction System



PCR System

Supplies and Reagents

Some materials may be usable on more than one system



Primers and Probes



Master Mix



Collection and Media



General Reagents

[Return to the Main Menu](#)


[Return to Overview of System Types](#)

[Additional Resources](#)

[Return to FDA FAQ](#)

PCR Testing Process

Schematic of SARS-CoV-2 Molecular Diagnostic Assay




Select the STEPs above for more information about the PCR testing process

[Return to the Main Menu](#)

[Additional Resources](#)

[Return to FDA FAQ](#)



INPUT FOR THIS STEP? A sample from the patient to test for the presence of SARS-CoV-2

WHAT SUPPLIES ARE NEEDED?

A swab to collect the sample.

Transport Media in a collection tube to transport the sample without degradation. The swab carrying the sample is placed in a collection tube filled with transport media.

WHAT IS THE OUTCOME? A sample for PCR analysis

[Return to the Main Menu](#)

[Return to PCR Test Process Overview](#)

[Additional Resources](#)

[Return to FDA FAQ](#)



INPUT FOR THIS STEP? A sample for PCR analysis

WHAT SUPPLIES ARE NEEDED?

Extraction System

Extraction Reagents

- Lysis Buffer
- RNA Extraction Control
- Human Specimen Control
- Positive Control specific to SARS-CoV-2


WHAT IS THE OUTCOME? Extracted nucleic acids

[Return to the Main Menu](#)

[Return to PCR Test Process Overview](#)

[Additional Resources](#)

[Return to FDA FAQ](#)



INPUT FOR THIS STEP? Extracted nucleic acid

WHAT SUPPLIES ARE NEEDED?

PCR REAGENTS

- PCR buffer
- dNTPs (building blocks of nucleic acids)
- Reverse Transcriptase (RT)
- No Template Control
- Polymerase (Enzyme)
- Positive Control (specific to SARS-CoV-2)
- Primers and Probe (SARS-CoV-2 specific)


WHAT IS THE OUTCOME? Fluorescence output signal indicating presence of SARS-CoV-2 RNA

[Return to the Main Menu](#)

[Return to PCR Test Process Overview](#)

[Additional Resources](#)

[Return to FDA FAQ](#)





INPUT FOR THIS STEP?

Fluorescence output signal

WHAT IS THE OUTCOME?

Report of output indicating presence of SARS-CoV-2 virus



[Return to the Main Menu](#)

[Return to PCR Test Process Overview](#)

[Additional Resources](#)

[Return to FDA FAQ](#)

SELECT A PCR TEST COMPONENT

Primers and Probes

Master Mix

PCR System

Controls

Extraction System/Kits

General Consumables

Substitution Options*

** Substitution Options are illustrated with the CDC test for this tool. Labs can identify and validate alternative options to other authorized assays.*

Select a Test Component on the left to explore possible substitutions for the Open Mix and Match style CDC molecular test

Return to the Main Menu

Additional Resources

Return to FDA FAQ

SELECT A PCR TEST COMPONENT

Primers and Probes



Master Mix

PCR System

Controls

Extraction System/Kits

General Consumables

Primers and Probes are the key reagents for direct detection and identification of a given target sequence.

Targets listed in the CDC EUA

➤ CDC N1

➤ CDC N2

➤ RNaseP

CDC Vetted Acceptable Manufacturers (per [CDC guidance](#))

- ✓ BioSearch Technologies
- ✓ Integrated DNA Technologies

Return to the Main Menu

Return to Substitution Options

Additional Resources

Return to FDA FAQ

SELECT A PCR TEST COMPONENT

Primers and Probes

Master Mix

PCR System

Controls

Extraction System/Kits

General Consumables

The master mix contains the reagents for the PCR process. Alternatives will need to be validated with the PCR system.

Master mixes listed in the CDC EUA

- Quantabio qScript XLT
- Quantabio UltraPlex
- Promega GoTaq
- Thermofisher TaqPath

Validation Needed

CAUTION

Equivalent Master Mix

CAUTION

Return to the Main Menu

Return to Substitution Options

Additional Resources

Return to FDA FAQ

SELECT A PCR TEST COMPONENT

Primers and Probes

Master Mix

PCR System

Controls

Extraction System/Kits

General Consumables

The PCR system is the instrument performing amplification.

PCR system listed in the CDC EUA

- Applied Biosystems 7500 Fast Dx Real-Time PCR

Alternatives (see [FDA FAQ*](#))

- ✓ Applied Biosystems 7500 Fast Real-Time PCR
- ✓ Applied Biosystems QuantStudio 6 Flex
- ✓ Applied Biosystems QuantStudio Dx
- ✓ QIAGEN Rotor-Gene Q MDx

** Check [FDA FAQ](#) for updates after 6/3/2020*

Return to the Main Menu

Return to Substitution Options

Additional Resources

Return to FDA FAQ

SELECT A PCR TEST COMPONENT

Primers and Probes

Master Mix

PCR System

Controls

Extraction System/Kits

General Consumables

Controls are the reagents that monitor the extraction and amplification process.

Controls listed in the CDC EUA

- CDC Positive Control
- CDC Extraction Control

Alternatives (see [FDA FAQ*](#))

- ✓ BEI Resources (N1, N2 Positive Control)
- ✓ Integrated DNA Technologies (N1, N2, Rnase P Positive Control)
- ✓ Human RNA (Rnase P Positive Control, Extraction Control)
- ✓ Twist Bioscience (Synthetic RNA Controls)

** Check [FDA FAQ](#) for updates after 6/3/2020*

Return to the Main Menu

Return to Substitution Options

Additional Resources

Return to FDA FAQ

SELECT A PCR TEST COMPONENT

Primers and Probes

Master Mix

PCR System

Controls

Extraction System/Kits

General Consumables

The Extraction system is the instrument performing extraction. The Extraction kit contains the reagents for the Extraction process.

Extraction Systems/kits listed in the CDC EUA

- QIAGEN/QIAamp DSP Viral RNA Mini Kit
- QIAGEN/QIAamp Viral RNA Mini Kit
- QIAGEN EZ1 Advanced XL/EZ1 DSP Virus Kit, Buffer AVL
- QIAGEN EZ1 Advanced XL/EZ1 Virus Mini Kit v2.0, Buffer AVL

Alternatives (see [FDA FAQ*](#))

- | | |
|---|--|
| ✓ QIAGEN QIAcube/QIAamp DSP Viral RNA Mini Kit | ✓ bioMérieux EMAG/EasyMAG Extraction Reagents, EMAG 1000µL Tips |
| ✓ QIAGEN QIAcube/ QIAamp Viral RNA Mini Kit | ✓ KingFisher Flex/Omega Bio-Tek Mag-Bind Viral DNA/RNA 96 Kit |
| ✓ Roche MagNA Pure LC/Total Nucleic Acid Kit | ✓ Applied Biosciences MagMAX Express/Viral/Pathogen Ultra Nucleic Acid Isolation Kit |
| ✓ Roche MagNA Pure Compact/Nucleic Acid Isolation Kit I | ✓ Promega Maxwell RSC 48/Maxwell RSC Viral Total Nucleic Acid Purification Kit |
| ✓ Roche MagNA Pure 96/DNA and Viral RNA Small Volume Kit | |
| ✓ bioMérieux NucliSENS easyMAG/EasyMAG Extraction Reagents, Biohit Pipette Tips | |
| ✓ Beckman RNAdvance Viral XP extraction kit | |

*** Check [FDA FAQ](#) for updates after 6/3/2020**

Return to the Main Menu

Return to Substitution Options

Additional Resources

Return to FDA FAQ

SELECT A PCR TEST COMPONENT

Primers and Probes

Master Mix

PCR System

Controls

Extraction System/Kits

General Consumables

General consumable laboratory supplies used in PCR

General consumables listed in the CDC EUA

- Molecular grade water, nuclease-free
- Aerosol barrier pipette tips
- Microcentrifuge tubes
- DNA Decontamination Reagent (DNAZap)
- Surface Decontaminant (RNAse Away)
- PCR Reaction Plates
- Vortex Mixer
- Microcentrifuge
- Micropipettes
- 10% Bleach
- MicroAmp Optical 8-cap Strips

Return to the Main Menu

Return to Substitution Options

Additional Resources

Return to FDA FAQ

Additional Resources

- [FDA Frequently Asked Questions on Diagnostic Testing for SARS-CoV-2](#)
- [CDC 2019-Novel Coronavirus \(2019-nCoV\) Real-Time RT-PCR Diagnostic Panel](#)
- [Enforcement Policy for Viral Transport Media During the Coronavirus Disease 2019 \(COVID-19\) Public Health Emergency - Guidance for Commercial Manufacturers, Clinical Laboratories, and Food and Drug Administration Staff \(July 2020\)](#)